


Elite Imaging Table



Precision 0047-TM5(H)

*Unobstructed 3-D
c-arm imaging*

500lbs/227kg weight capacity

High Quality Hospital Grade Mobile C-Arm Imaging Table

The "Elite" series of imaging tables are designed for high patient throughput utilizing the highest quality components and designed for ease-of-use by the medical staff with patient comfort in mind.

The tabletop has been designed to give maximum unobstructed imaging allowing 360° rotation radiation-free coverage around the longitudinal axis of the carbon fiber tabletop which makes it ideal for 3D imaging studies.

The tabletop pad is also designed to give the patient maximum comfort being manufactured utilizing our patient form-fitting Visco/Memory "Sofcare" High Density foam.

150mm (58") Full 360° Unencumbered Rotational Imaging Area



Model: 0047-TM3

Multi-purpose, 3-movement
Imaging Table

- Motorized Elevation
- $\pm 15^\circ$ Trendelenburg Tilt
- $\pm 15^\circ$ Lateral Tilt

Model: 0047-TM4

Extend-Retract, 5-movement
Imaging Table

- Motorized Elevation
- 15.75" (400mm) Motorized Longitudinal Travel
- $\pm 15^\circ$ Trendelenburg Tilt
- $\pm 15^\circ$ Lateral Tilt





Movement One

- Motorized Elevating Lift
- Lift Travel: (11.75") 300mm
- Optional Travel: (15.75") 400mm



Movement Two

- Motorized Trendelenburg Tilt
- ± 15 Degree Tilt

Movement Three

- Motorized Extend and Retract
- Longitudinal Travel: (15.75") 400mm



Movement Four

- Motorized Lateral Tilt
- Lateral Tilt: ± 15 Degree Tilt



Movement Five

- Motorized Lateral Drive
- Lateral Travel: ($\pm 4"$) ± 100 mm

Accessories



Wireless Hand Control



Foot Control for all Movements



Accessory Clamp



Patient Arm Board and Pad

ARCOMA-MEDSTONE MOBILE C-ARM IMAGING TABLES

0047-SERIES SPECIFICATIONS

		TM-5ER(H)	TM-4ER(H)	TM-4ER	TM-3ER	TM-3MP	TM-2ER	TL-2/TT-2	TE-1
(1) Motorized Elevation (Lift) Travel	Minimum	26" (66.5cm)	26" (66.5cm)	30" (76cm)	26" (66.5cm)	27" (66.5cm)	26" (69cm)	25" (64cm)	22" (56cm)
	Maximum	38" (96cm)	38" (96cm)	41.75" (106cm)	38" (96cm)	98cm (39")	42.75(109cm)	41.75 (106cm)	37.75" (96cm)
(2) Motorized Trendelenburg Tilt		± 15°	± 15°	± 15°	± 15°	NO	NO	± 15°	NO
(3) Motorized Longitudinal Table-top Travel		15¾" 40cm	15¾" 40cm	15¾" 40cm	NO	NO	15¾" 40cm	NO	NO
(4) Motorized Lateral Table-top Tilt		± 15°	± 15°	± 15°	± 15°	± 15°	NO	NO	NO
(5) Motorized Lateral Table-top Travel		± 4" 100mm	NO	NO	NO	NO	NO	NO	NO
(6) Table Weights		510lbs (232kg)	490lbs (223kg)	510lbs (232kg)	510lbs (232kg)	450lbs (205kg)	500lbs (227kg)	420lbs (191kg)	410lbs (186kg)
(5) Shipping Weights		700lbs (318kg)	680lbs (309kg)	700lbs (318kg)	700lbs (318kg)	640lbs (291kg)	690lbs (314kg)	610lbs (277kg)	600lbs (273kg)
(7) Shipping Crate Size, Imperial	88" x 31" x 41"								
(8) Shipping Crate Size, Metric	230cm x 79cm x 104cm								
Carbon Fiber Table-top 22" wide x 22" long, with Facial access cutout		Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Carbon Fiber Rectanaglar Table-top 22" wide x 82" long, with or without Facial access cutout		OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Patient Weight Capacity		500lbs 227kg	500lbs 227kg	500lbs 227kg	500lbs 227kg	500lbs 227kg	500lbs 227kg	500lbs 227kg	500lbs 227kg
Un-obstructed Table-top Imaging Area		58" 1473 mm	58" 1473 mm	58" 1473 mm	58" 1473 mm	58" 1473 mm	58" 1473 mm	58" 1473 mm	58" 1473 mm
Radiolucent Sofcare Pad		YES	YES	YES	YES	YES	YES	YES	YES
Pateint Restraints (3)		YES	YES	YES	YES	YES	YES	YES	YES

OPTIONS

Dual Patient Arm Boards with Pads	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Snap-on Arm Board with Mounting Clamp	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Wireless Hand Control with Receiver	NO	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
Foot Pedel Control	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL



Technical specifications are subject to change without notice.
© 2009 Arcoma North America, Inc. 0047-005-A (03-09)



North American Office 23112 Alcalde Dr., Suite A, Laguna Hills, CA 92653
Worldwide Headquarters Anavagen I, Vaxjo SE-325 46, Sweden